Thank you for purchasing IALUADAPTER®, which enables catheter-free instillation of the bladder. Compared to the traditional way of bladder instillation – which involves a catheter – using IALUADAPTER® has several advantages. The process is minimally invasive and entirely painless. A catheter, even if it is used in the most cautious way, can cause microlesions and be the source of catheter-associated infections. These issues are not experienced if IALUADAPTER® is applied. Moreover, IALUADAPTER® enables the treatment of the bladder and the urethra at the same time, which is challenging or impossible with a catheter.

#### Intended Use

Performing bladder instillation with any solution prescribed by a therapist without a catheter, in a safe, painless, and minimally invasive way. Generally used solutions include substances for regenerating the inner mucosal surface of the bladder - the GAG-layer - (e.g., hyaluronic acid, chondroitin sulphate, heparin, DMSO, or any combination of them), local anesthetics (e.g. lidocaine), antibiotics, antiseptics, certain hormones, chemotherapeutic agents (mitomycin C, epirubicin, gemcitabine, valrubicin, docetaxel) or immunotherapeutic agents (BCG - Bacillus Calmette-Guérin) to treat bladder cancer, and agents to treat hemorrhagic cystitis.

#### Indications

Indications for which IALUADAPTER® can be successfully applied are Interstitial Cystitis/Bladder Pain Syndrome (IC/BPS), recurrent Urinary Tract Infections (rUTIs), certain forms of Hemorrhagic Cystitis (Chemotherapy Induced Cystitis, Radiation Cystitis), Severe Overactive Bladder Syndrome (OAB), and local chemo- or immunotherapy of Bladder Cancer. Product could be used also with physiological saline, phase solutions, anti-inflammatory agents (NSAIDS, corticosteroids), beta-3 adrenergic agonists (mirabegron), anticholinergic agents (oxybutynin, tolterodine, solifenacin, trospium, darifenacin), and contrast agents for imaging. IALUADAPTER® can be applied to male or female patients. Product is intended for transient use.

#### Contraindications

Use of product is not advised on patients with known or suspected hypersensitivity to any of the materials employed.

#### Intended Users

IALUADAPTER® can be used by medical professionals who are expert in the field of bladder instillation and experienced in the use of the device. IALUADAPTER® can be used by any non-professional person providing they have been trained on how to use the device by a professional priority.

#### Materials Used

The product is made of medical grade thermoplastic elastomer (TPE).

#### Instructions for Use

1. Perform the instillation only in places with suitable hygiene level – poor hygienic environment may result in complications. Before starting the treatment, disinfect not only the urethral orifice but also the 3 cm surrounding area and the inner surface of the labia with any mucosal antiseptic solution.
2. Discard the product if the packaging is damaged. Unpack IALUADAPTER® and attach its connective tail (4) to the syringe. After taking the product out of the packaging, inspect the product visually for mechanical damage – do not use the product if there is any present. After unpacking, use the device immediately. IALUADAPTER® tail (4) can be attached to both Luer Slip and Luer Lock syringes easily. Although IALUADAPTER® is sterile, it might be worth applying disinfectant agent on it, from the tip (1) to the sealing collar (2). This process can prevent the drift of the bacteria into the bladder, and inactivate the bacteria located in the first couple of millimeters of the urethra.
3. Always make it sure that IALUADAPTER® is inserted into the urethra, and not into the vagina in case of female patients.

4. For exposing the orifice of female patients, the fingers should lie as close to the orifice as they can be. The exosition has to happen by pulling the labia apart laterally, too (and not only upwards).
5. The optimal instillation can be performed if the axes of the urethra and the syringe are aligned in parallel. In this case, the resistance can be experienced by using the syringe decreases characteristically.
6. Only the tip of IALUADAPTER® should penetrate the urethra. The isolating collar (2) has to cover the orifice properly and be pushed against the nearby structures gently.
7. Start performing the instillation. Given that the isolating collar (2) is covering the orifice properly, the instilled solution raises the intra-urethral pressure, which opens the sphincter, thus the solution enters the bladder.
8. If the sphincter becomes constricted due to pain or stress, tell the patient to relax. Wait until the constrictions end; the speed of the instillation should be decreased or paused.
9. In case of leakage, a part of the solution does not reach the urethra, it flows away instead. The two most frequent situations in which leakage happens are applying too low or too high pressure with IALUADAPTER®. Adjust the pressure accordingly; the ideal amount can be experienced as the sudden lowering in resistance of the plunger of the syringe.
10. In most cases, the leaking solution causes neither infections nor complications by entering the vagina or contacting the nearby structures. In case of doubt, consult the Instructions for Use or general description of the solution being instilled.

#### Precautions and Warnings

- IALUADAPTER® is sterile; it is for single use only. Only a professional person can prescribe the use of IALUADAPTER® for self-instillation. Use is allowed to professionals or priority trained individuals. Before starting a self-therapy with IALUADAPTER® always consult with your therapist about
- the way how the self-treatment should be performed
  - the solution you have to instill into your bladder
  - the frequency and the total number of self-treatments.

IALUADAPTER® is designed for administering solutions into the bladder, through the orifice of the urethra. Never use IALUADAPTER® for delivering any liquid into any other organ. IALUADAPTER® is intended for instilling solution into the bladder. Do not use it for draining of urine or any other liquid. Do not reuse or resterilize IALUADAPTER®.

#### Danger!

In case of reuse, risk of transmitting infections and occurrence of inflammation is present. Do not use the product if its packaging is damaged. After use, product has to be discarded in accordance with medical practice and local regulations. Keep IALUADAPTER® away from children younger than 3 years, since it is a suffocation hazard for them. In case of performing an instillation with the help of a non-professional person, ensure that they have read this Instructions for Use before starting the treatment, and have understood the instructions of your therapist fully. Ensure that you will not have to change posture during the self-treatment against your will; do not commit public indecency. In case of any side effects caused by the treatment or the instilled solution, contact your therapist immediately. In case of adverse effects caused directly by IALUADAPTER®, inform TIK d.o.o. and National Competent Authority on the incident. In case the product is not used in accordance with the instructions provided by TIK d.o.o., the manufacturer or the distributor do not take any responsibility for any injuries or other damage that may arise. For any complications or adverse effects caused by the solution instilled with IALUADAPTER® neither the manufacturer, nor the distributor take any responsibility.

#### Storage

If the product or its sterile packaging is damaged, the risk of transferring infections may be increased/introduced. Keep dry and store within the prescribed temperature range. Inappropriate storage can disrupt product sterility, affecting the device's safety and usability. The product should be stored in its original packaging in a dry place, in a temperature range from 5 to 40 °C (from 41 to 104 °F). The product should be protected from extensive heat or direct sunlight. The sterility for a properly stored product is guaranteed for 5 years unless the packaging is opened or damaged.

For further information please visit websites: [www.urosystem.com](http://www.urosystem.com), [www.urodapter.com](http://www.urodapter.com) and <https://app.urodapter.com>, where you can find a detailed Tips & Tricks guide of using IALUADAPTER® for therapists.

This Instructions for Use is a subject of updates. The latest version of the document is always available on the [www.urosystem.com](http://www.urosystem.com) webpage.

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Qarad Suisse S.A.