



Knowledge for life
since 1951.

Catalog




Knowledge for life since 1951.



**Knowledge for life
since 1951.**



Watch the video.




Select an European manufacturer and retain the flexibility of deliveries. Participate in the development of products for your needs.

TIK d.o.o. company is the right choice, if you are looking for a reliable supplier of high quality and user-friendly products for human medicine.

As an European manufacturer of medical devices with more than 70 years of tradition, we have always been faithful to the quality at all levels of business. Based on years of experience and our own development, we offer products according to the modern needs of the market. On account of our flexibility we can fulfill your needs and desires. Our products are following all quality related standards to help people in facilitating a healthy life.

You can be carefree. Cooperation with TIK d.o.o. brings you new reliability, professional support and inspires confidence in end users.



INDEX

7

UROLOGY

- 8 greencath® brand
- 10 greencath® - Coated Nelaton urethral catheter
- 11 greencath® PLUS - Coated Nelaton urethral catheter with water sachet
- 12 greencath® SOFT - Coated soft tip urethral catheter
- 13 greencath® SOFT PLUS - Coated soft tip urethral catheter with water sachet
- 14 greencath® Tiemann PLUS - Coated Tiemann urethral catheter with water sachet
- 14 greencath® Tiemann SOFT PLUS - Coated soft tip Tiemann urethral catheter with water sachet
- 15 greencath® BASIC - Coated Nelaton urethral catheter with water sachet

- 16 DO.TIK - Nelaton urethral catheter with polished eyes
- 16 DO.TIK SOFT - Soft tip urethral catheter with polished eyes

- 17 Nelaton urethral catheter
- 17 Nelaton LL urethral catheter

- 18 Tiemann urethral catheter
- 18 Connector between catheter and urine bag

- 19 Urine bag
- 19 Urine bag LIFE

21

RESPIRATORY THERAPY

- 22 Aspiration catheter
- 22 Aspiration catheter with vacuum control
- 23 Aspiration catheter with vacuum control - TIK PLUS
- 23 Aspiration catheter with curved tip and vacuum control

GASTROENTEROLOGY



- 26 Feeding tube with ENFit® connector
- 26 Feeding tube
- 27 Stomach tube
- 28 Duodenal tube
- 28 Rectal tube
- 29 Catheter stopper

25

GYNAECOLOGY

- 32 Follicle aspiration needle

31






UROLOGY

- **greencath**[®] - Coated Nelaton urethral catheter
- **greencath**[®] PLUS - Coated Nelaton urethral catheter with water sachet
- **greencath**[®] SOFT - Coated soft tip urethral catheter
- **greencath**[®] SOFT PLUS - Coated soft tip urethral catheter with water sachet
- **greencath**[®] Tiemann PLUS - Coated Tiemann urethral catheter with water sachet
- **greencath**[®] Tiemann SOFT PLUS - Coated soft tip Tiemann urethral catheter with water sachet
- **greencath**[®] BASIC - Coated Nelaton urethral catheter with water sachet
- DO.TIK - Nelaton urethral catheter with polished eyes
- DO.TIK SOFT - Soft tip urethral catheter with polished eyes
- Nelaton urethral catheter
- Nelaton LL urethral catheter
- Tiemann urethral catheter
- Connector between catheter and urine bag
- Urine bag
- Urine bag LIFE




Urethral catheters

greencath[®]



[Watch the video](#)

greencath[®] catheters are innovative, user friendly catheters, developed with an advanced technological process and with a special formula of TIK's hydrophilic coating. This special hydrophilic coating, which is a product of continuous development, reduces the friction between the catheter and the urethral mucosal lining, thus reducing stress on the organism.



SLEEVE

Some dimensions have an additional sleeve allowing aseptic handling.



SOFT TIP

Soft, flexible tip enables easy catheter introduction in case of enlarged prostate and various urethral masses.



HYDROPHILIC COATING

Catheters feature a biocompatible coating that is activated by water in 30 seconds, a longer drying time after coating activation, and their excellent sliding capacity. Its compactness and good adhesion to the base catheter ensures that the coating does not detach after activation with water.



POLISHED SIDE EYES

Polished side eyes significantly reduce the possibility of urethral injuries and consequential urinary tract infections.



TIEMANN TIP

Tiemann tip enables easy catheter introduction in case of enlarged prostate and various urethral masses.



WATER SACHET

Water sachet attached, for easy activation of the coating, containing specially purified, sterile water.



PRIMARY PACKAGING

The opening system is simple and practical. The catheter in its primary packaging can be fitted in the upright position with two attachment systems - using a two-sided sticker on the back of the packaging, or by hanging using a hole on the blister.



greencath®

Coated Nelaton urethral catheter

- For daily catheterisation, for short-term use for introduction through the urethra for urine drainage or therapeutic purposes.
- Hydrophilic coating reduces the friction between the catheter and urethral mucosal lining.
- Minimal residue of the hydrophilic coating in the urethra reduces stress to the organism caused by absorption of the hydrophilic coating after the use of the catheter.
- Polished side eyes significantly reduce the possibility of urethral injuries and consequential urinary tract infections.
- Connectors are colour coded for easier dimension identification.



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)	REMARK		
59140	■ blue	8	400	2,6	Male	30	600
59141	■ black	10	400	3,3	Male	30	600
59142	□ white	12	400	4,0	Male	30	600
59143	■ green	14	400	4,6	Male	30	600
59144	■ orange	16	400	5,3	Male	30	600
59145	■ red	18	400	6,0	Male	30	600
59146	■ blue	8	180	2,6	Female	30	900
59147	■ black	10	180	3,3	Female	30	900
59148	□ white	12	180	4,0	Female	30	900
59149	■ green	14	180	4,6	Female	30	900
59150	■ orange	16	180	5,3	Female	30	900
59151	■ red	18	180	6,0	Female	30	900
59152	■ green	6	280	2,0	Paediatric	30	600
59153	■ blue	8	280	2,6	Paediatric	30	600
59154	■ black	10	280	3,3	Paediatric	30	600

UROLOGY



Hydrophilic coating



Polished side eyes

greencath® PLUS

Coated Nelaton urethral catheter with water sachet

- For daily catheterisation, for short-term use for introduction through the urethra for urine drainage or therapeutic purposes.
- Hydrophilic coating reduces the friction between the catheter and urethral mucosal lining.
- Minimal residue of the hydrophilic coating in the urethra reduces stress to the organism caused by absorption of the hydrophilic coating after the use of the catheter.
- Polished side eyes and tube flexibility additionally reduce friction between the catheter and urethral mucosal lining.
- Water sachet attached, for easy activation of the coating, containing specially purified, sterile water.
- Some dimensions have an additional sleeve allowing aseptic handling.
- Connectors are colour coded for easier dimension identification.

STERILE EO



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)	REMARK		
62294*	■ blue	8	400	2,6	Male	30	600
62295*	■ black	10	400	3,3	Male	30	600
62296*	□ white	12	400	4,0	Male	30	600
62297*	■ green	14	400	4,6	Male	30	600
62298*	■ orange	16	400	5,3	Male	30	600
62299*	■ red	18	400	6,0	Male	30	600
62300	■ blue	8	180	2,6	Female	30	900
62301	■ black	10	180	3,3	Female	30	900
62302	□ white	12	180	4,0	Female	30	900
62303	■ green	14	180	4,6	Female	30	900
62304	■ orange	16	180	5,3	Female	30	900
62305	■ red	18	180	6,0	Female	30	900
62306*	■ green	6	280	2,0	Paediatric	30	600
62307*	■ blue	8	280	2,6	Paediatric	30	600
62308*	■ black	10	280	3,3	Paediatric	30	600

*with sleeve



Hydrophilic coating



Polished side eyes



Water sachet



Sleeve



greencath® SOFT

Coated soft tip urethral catheter

- For daily catheterisation, for short-term use for introduction through the urethra for urine drainage or therapeutic purposes.
- Hydrophilic coating reduces the friction between the catheter and urethral mucosal lining.
- Minimal residue of the hydrophilic coating in the urethra reduces stress to the organism caused by absorption of the hydrophilic coating after the use of the catheter.
- Polished side eyes significantly reduce the possibility of urethral injuries and consequential urinary tract infections.
- Soft, flexible tip enables easy catheter introduction in case of enlarged prostate and various urethral masses.
- Connectors are color coded easier dimension identification



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)	REMARK		
62549	■ blue	8	400	2,6	Male	30	600
62550	■ black	10	400	3,3	Male	30	600
62551	□ white	12	400	4,0	Male	30	600
62552	■ green	14	400	4,6	Male	30	600
62553	■ orange	16	400	5,3	Male	30	600
62554	■ blue	8	180	2,6	Female	30	900
62555	■ black	10	180	3,3	Female	30	900
62556	□ white	12	180	4,0	Female	30	900
62557	■ green	14	180	4,6	Female	30	900
62558	■ orange	16	180	5,3	Female	30	900
62560	■ blue	8	310	2,6	Paediatric	30	600
62561	■ black	10	310	3,3	Paediatric	30	600

UROLOGY



Hydrophilic coating



Polished side eyes



Soft tip

greencath® SOFT PLUS

Coated soft tip urethral catheter with water sachet

- For daily catheterisation, for short-term use for introduction through the urethra for urine drainage or therapeutic purposes.
- Hydrophilic coating reduces the friction between the catheter and urethral mucosal lining.
- Minimal residue of the hydrophilic coating in the urethra reduces stress to the organism caused by absorption of the hydrophilic coating after the use of the catheter.
- Polished side eyes significantly reduce the possibility of urethral injuries and consequential urinary tract infections.
- Soft, flexible tip enables easy catheter introduction in case of enlarged prostate and various urethral masses.
- Water sachet attached, for easy activation of the coating, containing specially purified, sterile water.
- Some dimensions have an additional sleeve allowing aseptic handling.



STERILE EO



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)	REMARK		
62562	■ blue	8	400	2,6	Male	30	600
62563*	■ black	10	400	3,3	Male	30	600
62564*	□ white	12	400	4,0	Male	30	600
62565*	■ green	14	400	4,6	Male	30	600
62566*	■ orange	16	400	5,3	Male	30	600
62567	■ blue	8	180	2,6	Female	30	900
62568	■ black	10	180	3,3	Female	30	900
62569	□ white	12	180	4,0	Female	30	900
62570	■ green	14	180	4,6	Female	30	900
62571	■ orange	16	180	5,3	Female	30	900
62573	■ blue	8	310	2,6	Paediatric	30	600
62574*	■ black	10	310	3,3	Paediatric	30	600

*with sleeve



Hydrophilic coating



Polished side eyes



Water sachet



Soft tip



Sleeve



greencath® Tiemann PLUS

Coated Tiemann urethral catheter with water sachet

- For daily catheterisation, for short-term use for introduction through the urethra for urine drainage or therapeutic purposes.
- Hydrophilic coating reduces the friction between the catheter and urethral mucosal lining.
- Minimal residue of the hydrophilic coating in the urethra reduces stress to the organism caused by absorption of the hydrophilic coating after the use of the catheter.
- Tiemann tip enables easy catheter introduction in case of enlarged prostate and various urethral masses.
- Water sachet attached, for easy activation of the coating, containing specially purified, sterile water.
- All dimensions have an additional sleeve allowing aseptic handling.
- Connectors are colour coded for easier dimension identification.



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)		
64337	■ black	10	400	3,3	30	600
64338	□ white	12	400	4,0	30	600
64339	■ green	14	400	4,6	30	600
64713	■ orange	16	400	5,3	30	600
64714	■ red	18	400	6,0	30	600



Hydrophilic coating



Tiemann tip



Water sachet



Sleeve

greencath® Tiemann SOFT PLUS

Coated soft tip Tiemann urethral catheter with water sachet

- For daily catheterisation, for short-term use for introduction through the urethra for urine drainage or therapeutic purposes.
- Hydrophilic coating reduces the friction between the catheter and urethral mucosal lining.
- Minimal residue of the hydrophilic coating in the urethra reduces stress to the organism caused by absorption of the hydrophilic coating after the use of the catheter.
- Soft, flexible Tiemann tip enables easy catheter introduction in case of enlarged prostate and various urethral masses.
- Water sachet attached, for easy activation of the coating, containing specially purified, sterile water.
- All dimensions have an additional sleeve allowing aseptic handling.
- Connectors are colour coded for easier dimension identification.



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)		
64763	■ black	10	400	3,3	30	600
64764	□ white	12	400	4,0	30	600
64765	■ green	14	400	4,6	30	600
64766	■ orange	16	400	5,3	30	600



Hydrophilic coating



Tiemann tip



Water sachet



Sleeve

greencath® BASIC

Coated Nelaton urethral catheter with water sachet

- For daily catheterisation, for short-term use for introduction through the urethra for urine drainage or therapeutic purposes.
- Hydrophilic coating reduces the friction between the catheter and urethral mucosal lining.
- Minimal residue of the hydrophilic coating in the urethra reduces stress to the organism caused by absorption of the hydrophilic coating after the use of the catheter.
- Side eyes significantly reduce the possibility of urethral injuries and consequential urinary tract infections.
- Water sachet attached, for easy activation of the coating, containing specially purified, sterile water.
- Connectors are colour coded for easier dimension identification.



STERILE | EO



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)	REMARK		
64460	□ white	12	400	4,0	Male	30	600
64461	□ white	12	180	4,0	Female	30	900



Hydrophilic coating



Water sachet



DO.TIK

Nelaton urethral catheter with polished eyes

- Matt surface of the tube allows introduction with minimal use of lubricants.
- Tube flexibility and the shape of the side eyes enable easy catheter introduction into the urethra.
- Some dimensions have an additional sleeve allowing aseptic handling.
- Connectors are colour coded for easier dimension identification.



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)	REMARK		
61209	■ green	6	400	2,0	Male	60	840
61210	■ blue	8	400	2,6	Male	60	840
61211	■ black	10	400	3,3	Male	60	840
61212*	□ white	12	400	4,0	Male	60	840
61213*	■ green	14	400	4,6	Male	60	840
61214	■ orange	16	400	5,3	Male	60	840
61215	■ red	18	400	6,0	Male	60	840
61218	■ green	6	180	2,0	Female	60	1920
61219	■ blue	8	180	2,6	Female	60	1920
61220	■ black	10	180	3,3	Female	60	1920
61221	□ white	12	180	4,0	Female	60	1920
61222	■ green	14	180	4,6	Female	60	1920
61223	■ orange	16	180	5,3	Female	60	1920
61224	■ red	18	180	6,0	Female	60	1920

*with sleeve



Polished side eyes



Sleeve

DO.TIK SOFT

Soft tip urethral catheter with polished eyes

- Matt surface of the tube allows introduction with minimal use of lubricants.
- Tube flexibility and the shape of the side eyes enable easy catheter introduction into the urethra.
- Soft, flexible tip enables easy catheter introduction in case of enlarged prostate and various urethral masses.
- Suitable also for intermittent self-catheterization with Mitrofanoff bladder reconstruction.
- Some dimensions have an additional sleeve allowing aseptic handling
- Connectors are colour coded for easier dimension identification.



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)	REMARK		
62508	■ blue	8	400	2,6	Male	60	840
62509	■ black	10	400	3,3	Male	60	840
62510*	□ white	12	400	4,0	Male	60	840
62511*	■ green	14	400	4,6	Male	60	840
62512	■ orange	16	400	5,3	Male	60	840
62513	■ blue	8	180	2,6	Female	60	1920
62514	■ black	10	180	3,3	Female	60	1920
62515	□ white	12	180	4,0	Female	60	1920
62516	■ green	14	180	4,6	Female	60	1920
62517	■ orange	16	180	5,3	Female	60	1920
62519	■ blue	8	310	2,6	Paediatric	60	840
62520	■ black	10	310	3,3	Paediatric	60	840

*with sleeve



Polished side eyes



Soft tip



Sleeve

NELATON URETHRAL CATHETER

- For daily catheterisation, for short-term use for introduction through the urethra for urine drainage or therapeutic purposes.
- Matt surface of the tube allows introduction with minimal use of lubricants.
- Flexible tube and special shape of side eyes enable atraumatic insertion of the catheter into urethra.
- Shape of lateral eyes enables efficient drainage.
- Connectors are colour coded for easier dimension identification.

STERILE | EO



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)	REMARK		
54802	■ green	6	400	2,0	Male	60	840
54803	■ blue	8	400	2,6	Male	60	840
54804	■ black	10	400	3,3	Male	60	840
54805	□ white	12	400	4,0	Male	60	840
54806	■ green	14	400	4,6	Male	60	840
54807	■ orange	16	400	5,3	Male	60	840
54808	■ red	18	400	6,0	Male	60	840
54809	■ yellow	20	400	6,6	Male	60	840
54810	■ violet	22	400	7,3	Male	60	840
54820	■ green	6	180	2,0	Female	60	1920
54821	■ blue	8	180	2,6	Female	60	1920
54822	■ black	10	180	3,3	Female	60	1920
54823	□ white	12	180	4,0	Female	60	1920
54824	■ green	14	180	4,6	Female	60	1920
54825	■ orange	16	180	5,3	Female	60	1920
54826	■ red	18	180	6,0	Female	60	1920



NELATON LL URETHRAL CATHETER

- For taking urine sample with a syringe.
- Luer lock connector enables connecting with a syringe.
- Connectors are colour coded for easier dimension identification.

STERILE | EO



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)	REMARK		
53697	■ green	6	400	2,0	Male	60	840
54435	■ blue	8	400	2,6	Male	60	840
54436	■ black	10	400	3,3	Male	60	840





TIEMANN URETHRAL CATHETER

- For daily catheterisation, for short-term use for introduction through the urethra for urine drainage or therapeutic purposes.
- Matt surface of the tube allows introduction with minimal use of lubricants.
- Appropriate for users with enlarged prostate and various urethral masses.
- Connectors are colour coded for easier dimension identification.

STERILE EO



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)		
54812	■ blue	8	400	2,6	60	840
54813	■ black	10	400	3,3	60	840
54814	□ white	12	400	4,0	60	840
54815	■ green	14	400	4,6	60	840
54816	■ orange	16	400	5,3	60	840
54817	■ red	18	400	6,0	60	840
54818	■ yellow	20	400	6,6	60	840
54819	■ violet	22	400	7,3	60	840



Tiemann tip

CONNECTOR BETWEEN CATHETER AND URINE BAG

- Luer lock connector enables connection with a catheter which has Luer lock connector.
- Connector enables connection with urine bag.

STERILE EO



REF		
53688	60	660



URINE BAG

- Enables filling up to 800 mL.
- Luer lock connector can be used in combination with nephrostomy catheters.
- Adapter connector can be used in combination with indwelling urethral catheters.
- Integrated non-return valve and "T-release".
- Membrane for needlefree collection of urine sample.
- Plastic holder.
- Indicative scale in millilitres.
- Period of use up to 7 days.



REF	CONNECTOR TYPE	TUBE LENGTH (mm)	NOTE		
53889	Adapter	900	For indwelling urethral catheters	5	60
53900	Luer Lock	900	For nephrostomy catheters	5	60



URINE BAG LIFE

- Enables filling up to 800 mL.
- Luer lock connector can be used in combination with nephrostomy catheters.
- Adapter connector can be used in combination with indwelling urethral catheters.
- Integrated non-return valve and "T-release".
- Indicative scale in millilitres.
- Period of use up to 7 days.



REF	CONNECTOR TYPE	TUBE LENGTH (mm)	NOTE		
64454	Adapter	650	For indwelling urethral catheters	10	120
64455	Luer Lock	650	For nephrostomy catheters	10	120





RESPIRATORY THERAPY

- Aspiration catheter
- Aspiration catheter with vacuum control
- Aspiration catheter with vacuum control - TIK PLUS
- Aspiration catheter with curved tip and vacuum control



ASPIRATION CATHETER

- For suction from nose, throat or lungs.
- Catheter with distal eye and 2 side eyes.
- Centimeter marks.
- Connectors are colour coded for easier dimension identification.



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)		
54843	grey	5	300	1,7	60	840
54844	green	6	300	2,0	60	840
54845	blue	8	300	2,6	60	840
54846	black	10	300	3,3	60	840
64705	white	12	300	4,0	60	840
64706	green	14	300	4,6	60	840
64707	orange	16	300	5,3	60	840
54847	green	6	530	2,0	60	600
54848	blue	8	530	2,6	60	600
54849	black	10	530	3,3	60	600
54850	white	12	530	4,0	60	600
54851	green	14	530	4,6	60	600
54852	orange	16	530	5,3	60	600
54853	red	18	530	6,0	60	600
54854	yellow	20	530	6,6	60	600
56419	green	6	600	2,0	60	600
56420	blue	8	600	2,6	60	600
56421	black	10	600	3,3	60	600
56422	white	12	600	4,0	60	600
56423	green	14	600	4,6	60	600
56424	orange	16	600	5,3	60	600
56425	red	18	600	6,0	60	600



ASPIRATION CATHETER WITH VACUUM CONTROL

- For suction from nose, throat or lungs.
- Catheter with distal eye and 2 side eyes.
- Centimeter marks.
- Connector enables optimal vacuum control.
- Connectors are colour coded for easier dimension identification.



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)		
55126	grey	5	300	1,7	60	840
55127	green	6	300	2,0	60	840
55128	blue	8	300	2,6	60	840
55129	black	10	300	3,3	60	840
55114	grey	5	500	1,7	60	600
55115	green	6	500	2,0	60	600
55116	blue	8	500	2,6	60	600
55117	black	10	500	3,3	60	600
55118	white	12	500	4,0	60	600
55119	green	14	500	4,6	60	600
55120	orange	16	500	5,3	60	600
55121	red	18	500	6,0	60	600
56428	blue	8	550	2,6	60	600
56429	black	10	550	3,3	60	600
56430	white	12	550	4,0	60	600
56431	green	14	550	4,6	60	600
56432	orange	16	550	5,3	60	600
56433	red	18	550	6,0	60	600



ASPIRATION CATHETER WITH VACUUM CONTROL - TIK PLUS

- For suction from nose, throat or lungs.
- Smaller side eyes enable better suction through distal eye.
- Catheter with distal eye and 2 side eyes.
- Centimeter marks.
- Connector are colour coded for easier dimension identification.



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)		
62080	■ blue	8	530	2,6	60	600
62081	■ black	10	530	3,3	60	600
62082	□ white	12	530	4,0	60	600
62083	■ green	14	530	4,6	60	600
62084	■ orange	16	530	5,3	60	600
62085	■ red	18	530	6,0	60	600




ASPIRATION CATHETER WITH CURVED TIP AND VACUUM CONTROL

- For endobronchial suction procedures.
- Curved tip allows easy access to the left bronchus.
- 4 side eyes.
- Centimeter marks.
- Connectors are colour coded for easier dimension identification.



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)		
62473	■ black	10	530	3,3	60	600
62474	□ white	12	530	4,0	60	600
62475	■ green	14	530	4,6	60	600
62476	■ orange	16	530	5,3	60	600





GASTROENTEROLOGY

- Feeding tube with ENFit® connector
- Feeding tube
- Stomach tube
- Duodenal tube
- Rectal tube
- Catheter stopper




FEEDING TUBE WITH ENFit® CONNECTOR

- For enteral feeding in paediatrics.
- The tube only connects to syringes or systems that have a female ENFit® connector according to ISO 80369-3.
- 2 side eyes.
- Centimeter marks.
- Soft, rounded tip of the tube enables atraumatic insertion.
- Matted surface of the tube allows minimal use of lubricants.
- Connecting elements are colour-coded for quick and easy identification of CH.

STERILE EO



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)		
63873	grey	5	350	1,7	30	360
63874	green	6	350	2,0	30	360
63875	blue	8	350	2,6	30	360
63876	black	10	350	3,3	30	360
63877	grey	5	500	1,7	30	360
63878	green	6	500	2,0	30	360
63879	blue	8	500	2,6	30	360
63880	black	10	500	3,3	30	360
63881	white	12	500	4,0	30	360
63882	grey	5	1000	1,7	30	360
63883	green	6	1000	2,0	30	360
63884	blue	8	1000	2,6	30	360
63885	black	10	1000	3,3	30	360
63886	white	12	1000	4,0	30	360
64043	orange	16	1000	5,3	30	360
64044	red	18	1000	6,0	30	360



FEEDING TUBE

- For enteral feeding in paediatrics.
- Matted surface eliminates need for lubricant.
- Connector fits to luer syringes or administration lines (pumps).
- Centimeter marks.
- 2 side eyes.
- Feeding tube length 350 and 500 mm with 1 cm depth markings.
- Feeding tube length 1000 mm with 25 cm depth markings.
- Colour-coded connectors for easy identification of dimension.

STERILE EO



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)		
56164	grey	5	350	1,7	60	840
56165	green	6	350	2,0	60	840
56166	blue	8	350	2,6	60	840
56167	black	10	350	3,3	60	840
54861	grey	5	500	1,7	60	720
54862	green	6	500	2,0	60	720
54863	blue	8	500	2,6	60	720
54864	black	10	500	3,3	60	720
56104	white	12	500	4,0	60	720
56105	grey	5	1000	1,7	30	360
54865	green	6	1000	2,0	30	360
54866	blue	8	1000	2,6	30	360
54867	black	10	1000	3,3	30	360
56106	white	12	1000	4,0	30	360



STOMACH TUBE

- For feeding or emptying of stomach contents or gastric lavage.
- Colour-coded connectors for easy identification of dimension.
- Tube with closed rounded tip.
- 2 side eyes.
- Centimeter marks.
- Stomach tube length 800 mm with 10 cm depth markings: 45, 55 and 65 cm from the distal end.
- Stomach tube length 1200 mm with 10 cm depth markings: 45, 55, 65 and 75 cm from the distal end.

STERILE EO



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)		
54868	■ blue	8	800	2,6	25	300
54869	■ black	10	800	3,3	25	300
54870	□ white	12	800	4,0	25	300
54871	■ green	14	800	4,6	25	300
54872	■ orange	16	800	5,3	25	300
54873	■ red	18	800	6,0	25	300
54874	■ yellow	20	800	6,6	25	300
54875	■ violet	22	800	7,3	10	120
54876	■ blue	24	800	8,3	10	120
54877	□ white	26	800	8,6	10	120
54878	■ green	28	800	9,3	10	120
54879	■ grey	30	800	10,0	10	120
54880	■ black	35	800	11,6	10	120
54881	■ yellow	20	1200	6,6	10	120
54882	■ violet	22	1200	7,3	10	120





DUODENAL TUBE

- For aspiration of intestinal fluid in therapeutic or diagnostic purposes.
- Tube with closed tip.
- 4 side eyes.
- Colour-coded connectors for easy identification of dimension.
- Duodenal tube with 10 cm depth markings: 45, 55, 65 and 75 cm from the distal end.

STERILE EO



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)		
54855	■ blue	8	1200	2,6	25	300
54856	■ black	10	1200	3,3	25	300
54857	□ white	12	1200	4,0	25	300
54858	■ green	14	1200	4,6	25	300
54859	■ orange	16	1200	5,3	25	300
54860	■ red	18	1200	6,0	25	300



RECTAL TUBE

- For therapeutic purposes in case of blockade or cleaning of the bowel or for drainage of flatulence.
- 2 side eyes.
- Colour-coded connectors for easy identification of dimension.

STERILE EO



REF	COLOUR	CH	LENGTH (mm)	OUTER DIAMETER (mm)		
54827	■ red	18	400	6,0	50	700
54828	■ yellow	20	400	6,6	50	700
54829	■ violet	22	400	7,3	50	700
54830	■ blue	25	400	8,3	50	700
54831	■ green	28	400	9,3	50	700
54832	■ grey	30	400	10,0	50	700
54833	■ brown	32	400	10,6	50	700
54834	■ black	35	400	11,6	50	700



CATHETER STOPPER


- For closing of the connectors on catheters or tubes.

STERILE EO



REF		
2590	40	4000





GYNAECOLOGY

- Follicle aspiration needle



FOLLICLE ASPIRATION NEEDLE

- Follicle aspiration needle is used for ultrasonic guided, transvaginal harvesting of oocytes from follicles.
- Ergonomically designed needle holder.
- Atraumatic echo marker.
- Smooth inner surface of cannula and connecting tubes.
- Female Luer lock connector enables connecting to an aspirator or syringe.
- Transparent, kink resistant connecting tubes.
- Two bevel types available.
- Conical shaped tubing endings enable safe outflow into the tube.

STERILE EO



REF	OUTER DIAMETER OF THE NEEDLE (G/mm)	INNER DIAMETER OF THE NEEDLE (mm)	NEEDLE LENGTH (mm)	LENGTH OF TUBING FROM NEEDLE TO PLUG (mm)	LENGTH OF TUBING FROM PLUG TO ASPIRATOR (mm)	ECHO MARKER	ADDITIONAL NEEDLE FOR FLUSHING	BEVEL TYPE		
58831	17/1,4	1,18	350	670	700	✓	x	Lancet	10	120
59123	17/1,5	1,20	320	970	500	✓	x	Lancet	10	120
59436	17/1,5	1,20	320	670	700	✓	✓	Lancet	10	120
59587	17/1,5	1,20	240	670	700	✓	x	Lancet	10	120
59588	17/1,5	1,20	320	670	700	✓	x	Lancet	10	120
63393	17/1,4	0,92	320	970	1000	✓	x	Back cut bevel	10	120
63781	17/1,4	0,92	320	670	700	✓	x	Back cut bevel	10	120



GYNAECOLOGY



Watch the video



Explanation of the symbols that appear in the catalog:



Single use



Contains PVC



PVC free



Do not resterilize



Does not contain
phthalates
(DEHP-FREE)



Latex free



Sterile using
ethylene oxide



TIK d.o.o.

Goriška cesta 5b
5222 Kobarid,
Slovenia

T: +386 (0)5 389 07 00

@: info@tik.si

W: www.tik.si

